



Polylimb

*our sheets
mould consistently time
after time.*

Our knowledge of plastics has led to the development of materials made specifically for use in the orthotics and prosthetic profession.

Our business focus is on providing materials that can be relied to mould consistently and to offer top class technical backup and support giving our customers confidence in their manufacturing process and their products.



OUR RANGE INCLUDES-

- Ortho - AL
- Orthoreflex
- Orthoultraflex
- Orthoflex
- Orthoflex - 001

Polylimb Ortho-AI

WHAT MAKES US SPECIAL?

- Easy to fabricate, seam, drape and blister form
- Light weight
- Versatile
- Rigid
- High impact strength
- Good chemical resistance
- Good formability
- Maintain impact strength down to -200 C
- Has the ability to constantly flex without fracturing
- Shrinkage of less than 1%

APPLICATIONS

- Ankle Foot Orthosis - AFOs
- Bodyjackets
- Dynamic Ankle-Foot Orthosis DAFOs
- Definitive Prosthetic Sockets
- Functional Foot Orthosis - FFOs
- Knee Ankle Foot Orthosis - KAFOs
- Lightweight Prosthetic Liner Sockets

MATERIAL SPECIFICATIONS

Sheet size for orthotics :

1000 x 2000 mm

other sizes available on demand with MOQ

Available thickness :

2mm, 3mm, 4mm, 5mm, 6mm, 8mm

thickness upto 30mm available on demand with MOQ

Available Colours :

Natural

other colours available on demand with MOQ



Moulding temperature :

Orthotics : 180 C - 200 C

Prosthetic : 200 C - 220 C

Heating Time :

3 minutes per mm thickness

N.B. : As all ovens have different heating and sensing characteristics. Mentioned temperature and times of moulding is indicative only. Actual conditions will vary



| Polylimb Orthoreflex

WHAT MAKES US SPECIAL?

- Light weight
- More Rigid than Ortho - AL
- Good Toughness
- Good Impact strength
- Availability to flex without fracturing
- Shrinkage of less than 1%

APPLICATIONS

- Ankle Foot Orthosis - AFOs
- Dynamic Ankle-Foot Orthosis - DAFOs
- Definitive Prosthetic Sockets
- Frame Sockets
- Knee Ankle Foot Orthosis - KAFOs

MATERIAL SPECIFICATIONS

Sheet size for orthotics :

1000 x 2000 mm

other sizes available on demand with MOQ

Available thickness :

2mm, 3mm, 4mm, 5mm, 6mm, 8mm

thickness upto 30mm available on demand with MOQ

Available Colours :

Natural

other colours available on demand with MOQ



Moulding temperature :

Orthotics : 180 C - 200 C

Prosthetic : 200 C - 220 C

Heating Time :

3 minutes per mm thickness

N.B. : As all ovens have different heating and sensing characteristics. Mentioned temperature and times of moulding is indicative only. Actual conditions will vary



| Polylimb Orthoultraflex

WHAT MAKES US SPECIAL?

- Highly flexible thermoplastic sheet
- Excellent tear resistance, can be stitched or riveted
- No shrinkage
- The latest generation of polyolefin elastomer
- The lightest thermoplastic sheet material in the market



APPLICATIONS

- Flexible Liner Sockets
- Lightweight Prosthetic Liner Sockets
- Wrist splints
- Protective Helmets
- Neck Braces
- Body Jackets

MATERIAL SPECIFICATIONS

Sheet size for orthotics :

1000 x 2000 mm

other sizes available on demand with MOQ

Available thickness :

2mm, 3mm, 4mm, 5mm, 6mm, 8mm

thickness upto 30mm available on demand with MOQ

Available Colours :

Natural

other colours available on demand with MOQ

Moulding temperature :

Orthotics : 150 C - 180 C

Prosthetic : 150 C - 180 C

Heating Time :

3 minutes per mm thickness

N.B. : As all ovens have different heating and sensing characteristics. Mentioned temperature and times of moulding is indicative only. Actual conditions will vary recommended timings as a guide only.

| Polylimb Orthoflex-001

WHAT MAKES US SPECIAL?

- High molecular Weight Polyethylene
- Very Low Shrinkage
- Good chemical resistance
- Rigid
- Resists cracking
- Easy to fabricate

APPLICATIONS

- Ankle Foot Orthosis - AFOs
- Knee Ankle Foot Orthosis - KAFOs
- Functional Foot Orthosis - FFOs
- Body Jackets

MATERIAL SPECIFICATIONS

Sheet size for orthotics :

1000 x 2000 mm

other sizes available on demand with MOQ

Available thickness :

2mm, 3mm, 4mm, 5mm, 6mm, 8mm

thickness upto 30mm available on demand with MOQ

Available Colours :

Natural

other colours available on demand with MOQ



Moulding temperature :
Orthotics : 160 C - 180 C

Heating Time :
3 minutes per mm thickness

N.B. : As all ovens have different heating and sensing characteristics. Mentioned temperature and times of moulding is indicative only. Actual conditions will vary.



Polylimb Orthoflex

WHAT MAKES US SPECIAL?

APPLICATIONS

- Ankle Foot Orthosis - AFOs
- Body Jackets
- Wrist splints
- Neck braces
- Light weight AK/BK Sockets

MATERIAL SPECIFICATIONS

Sheet size for orthotics :

1000 x 2000 mm

other sizes available on demand with MOQ

Available thickness :

2mm, 3mm, 4mm, 5mm, 6mm, 8mm

thickness upto 30mm available on demand with MOQ

Available Colours :

Natural

other colours available on demand with MOQ



Moulding temperature :
Orthotics : 160 C - 180 C

Heating Time :
3 minutes per mm thickness

N.B. : As all ovens have different heating and sensing characteristics. Mentioned temperature and times of moulding is indicative only. Actual conditions will vary.

